


Biolog Microstation



[Link to Website](#)

Features:

Semi-Automated System for Microbial Identification
Turbidity Standards Included

Hardware:

Microstation ELx808BLG (S/N: 15050514)
Turbidimeter (S/N: 150367560)
Electronic Multichannel Pipette
Hewlett Packard ProDesk 400 G1MT (S/N: MXL4232213)
Hewlett Packard ProDisplay P17A (S/N: 3CQ450038K)
Brother HL-1110 Printer (S/N: E72141C5N556841)

Software:

MicroLog 3 MicroStation System Microbial Identification Software Ver. 5.2.2b35 (Pre-Installed & CD-ROM)

Retrospect Ver. 2.1.2 (Pre-Installed & CD-ROM)

Gen III Database Ver. 2.7.1 (Pre-Installed & CD-ROM)*

AN Database Ver. 6.01 (Pre-Installed & CD-ROM)*

Windows 7 Professional SP1 Operating System (Pre-Installed & DVD-ROM)

Driver Recovery DVD-ROM

HP ProDisplay P17A Monitor DVD-ROM

Brother Installer for Drivers Utilities and Users Manual CD-ROM

**Databases are single user agreement*

Manuals:

Electronic User Guide Pre-Installed

Ovation Multichannel Pipette Quick Reference Guide

Biolog Turbidimeter Owner's Manual

Additional Spare Parts and Consumables:

ITEM	MODEL/ PART No.	MANUFACTURER	SERIAL/ LOT No.	PACK SIZE	QUANTITY INCLUDED
25mL Reservoirs	3102	Biolog	25109	100	100
Innoculatorz, Individually Wrapped	3323	Biolog	4295/ 4152	100	200
96 Well Plate	-	Biolog	-	1	15
Well Plate Lid	-	Biolog	-	1	6
1250µL Pipette Tips, Sterile	4060-3132	VistaRak	15215P	192	1152

Also includes IQ/OQ/PQ documentation.